PHOSPHOLIN ES

Activated Partial Thromboplastin Time

English Version

INTENDED USE

Phospholin ES APTT reagent is an ellagic acid based reagent with soybean phospholipids, buffers, stabilizers and preservatives. Phospholin ES is intended for use as an activated partial thromboplastin time (APTT) reagent. The APTT test is a qualitative assay used in routine coagulation screening of patient plasma to detect deficiencies in the intrinsic pathway. It is also used to monitor heparin therapy and in the detection of Lupus Anticoagulants. Phospholin ES is to be used by qualified laboratory personnel.

SUMMARY AND PRINCIPLE

The activated partial thromboplastin time (APTT) is used to detect disorders in the intrinsic coagulation system, which involves coagulation factors II, V, VIII IX, XI, and XII. A modified version of the APTT is often used in assay systems that quantitate these factors. The APTT is commonly used for pre-surgical screening for intrinsic factor deficiencies (1), for monitoring heparin therapy (2), and in the detection of Lupus Anticoagulants (3).

In the basic screening test, the activated partial thromboplastin time indirectly measures the formation of thrombin by its action on fibrinogen resulting in a fibrin clot. In the test, citrated test plasma is mixed with APTT reagent for a specified period of time (typically 3 minutes) at 37°C followed by the addition of pre-warmed (37°C) calcium chloride (0.025 M). The clotting time is initiated by the addition of calcium chloride. The time (in seconds) required for clot formation is the activated partial thromboplastin time (APTT). Clot detection can be by mechanical, manual, or photo-optical measurement

REAGENTS

WARNING: FOR IN-VITRO DIAGNOSTIC USE ONLY
1. PHOSPHOLIN ES Reagent

Ingredients: The reagent contains 0.1mM ellagic acid with phospholipids derived from soybean lecithin. Buffers, stabilizers, and preservatives have been added.

Preparation for Use: The reagent is packaged ready for use. The vial should be inverted several times until a homogeneous suspension is obtained.

Storage and Stability: **PHOSPHOLIN ES** Reagent should be stored at 2 to 8°C when not in use and is stable until the date indicated on the vial.

DO NOT FREEZE.

Signs of Deterioration: The reagent is normally a homogeneous solution. A sediment may form upon standing that should dissipate easily upon mixing by inversion. Failure of normal plasma or controls to fall within established laboratory quality control ranges may indicate product deterioration.

2. Calcium Chloride Solution

Ingredients: 0.025M Calcium Chloride

Preparation for Use: The reagent is packaged ready for use.
Storage and Stability: 0.025M Calcium Chloride solution is stable until the date indicated on the vial.

Instruments

APTT times may be determined with **PHOSPHOLIN ES** using semi-automated and automated coagulation instruments.

Specimen Collection And Handling

NOTE: After initial whole blood collection, during testing all test tubes, syringes and pipettes should be plastic.

Specimen: Plasma obtained from whole blood anti-coagulated with 0.1 M sodium citrate. Use of hemolyzed, lipemic, or icteric samples should be avoided as these conditions may affect results, especially when using photo-optical instruments.

Specimen Collection: Nine parts of freshly collected whole blood should be immediately added to one part anticoagulant.

Specimen Preparation: Centrifuge the whole blood specimen at 2500 x g for 15 minutes (NCCLS H21-A2, 1991). Immediately separate the plasma from the red blood cells using a plastic pipette and place it in a plastic test tube. Perform the activated partial thromboplastin time assay within 2

hours.

Storage and Stability: Before and during testing the plasma sample should be maintained in the plastic test tubes at 2 to 8°C to insure further stability. If testing is delayed for more than 2 hours, plasma may be stored at -20°C or below for up to one month. Frozen samples should be thawed rapidly at 37°C before testing.

Materials Provided: Materials needed for **PHOSPHOLIN ES** assays are provided in the following packaging configurations.
PHOSPHOLIN ES

Reagent 5 x 5 mL vials, Phospholin ES
Reagent 5 x 5 mL vials, 0.025M Calcium Chloride
Reagent 5 x 10 mL vials, Phospholin ES
Reagent 5 x 10 mL vials, 0.025M Calcium Chloride

Materials and Equipment Required but not Provided:

Coagulation Instrument or 37°C water bath and timer
Reaction Cups or plastic test tubes
Pipettes to deliver 1.0 and 0.1 mL
Centrifuge
Distilled or deionized water
Control Plasmas:

PlasmaCon N
PlasmaCon L-1
PlasmaCon L-2

Test Procedure for APTT

NOTE: Throughout testing all test tubes, syringes and pi-pettes should be plastic.

I. Automated and Semi-Automated Methods
If using an instrument to perform this test, refer to the appropriate instrument Operator's Manual for instructions.

II. Manual method

- Collect blood specimen according to directions in SPECIMEN COLLECTION AND HANDLING SECTION.
- Centrifuge the anti-coagulated whole blood specimen at 2500 x g for 15 minutes or equivalent force-time.
- While the blood specimen is centrifuging, reconstitute the control plasma according to the package insert included with each control.
- Immediately after centrifugation, separate the plasma from the red blood cells and place in a plastic tube at 2 to 8°C until assayed. The maximum storage time at 2 to 8°C is 2 hours.
- Place the 0.025 M Calcium Chloride Reagent into a test tube and pre-warm to 37°C (requires approximately 5 minutes).
- Pipette 100 µL of the patient plasma or control plasma into a reaction tube.
- Gently mix the **PHOSPHOLIN ES** Reagent by inversion to re-suspend any sediment.
- Pipette 100 µL of **PHOSPHOLIN ES** Reagent into the reaction tube containing patient plasma or control plasma.
- Incubate the **PHOSPHOLIN ES** Reagent and patient plasma at 37°C for 3-5 MINUTES.
- Add 100 µL of pre-warmed 0.025 M Calcium Chloride, simultaneously starting a timer, and record the time (in seconds) required for clot formation.
- The test results (clot time) are reported directly in seconds as the APTT Time.

Heparin Monitoring

When using the APTT test to monitor heparin therapy, it is important to construct an *in-vitro* reference curve that reflects the average heparin response, since individual patients respond differently to heparin. In general, the accepted therapeutic range for heparin is 0.3 to 0.7 units/mL. The following precautions should be considered when monitoring heparin therapy.

- Time of collection is important, since heparin has an *in-vivo* half-life of only 1.5 hours.
- A base line APTT on each patient should be established before therapy starts in order to determine the degree of the patient APTT response. This should be correlated to the heparin dose response curve established by the testing

laboratory.

3. Heparin dose response curves should be constructed using the same heparin employed in therapeutic use to eliminate variables connected with heparin preparations derived from different sources (e.g., porcine mucosa or bovine lung).

4. Heparin dose response curves must be reestablished with each new lot number of reagent.

Quality Control

Each laboratory should establish a quality control program that includes normal and abnormal controls to evaluate instrument, reagent and technologist performance. The normal and abnormal controls should be tested daily prior to performing tests on patient plasmas. Monthly quality control charts (Levi Jennings) are recommended to determine the mean and standard deviation of control plasma.

A normal control such as PlasmaCon N and abnormal level 1 and level 2 controls such as PlasmaCon L-1 and PlasmaCon L-2 are recommended. If the controls do not perform within their reference ranges, patient results should be considered invalid and not reported.

Results

Results of the activated partial thromboplastin time testing should be reported as the APTT in seconds. These results should be interpreted in relation to the normal range for APTT testing in each laboratory. Times that are shorter or longer than the normal range may be indicative of an abnormal condition in the patients coagulation system.

Limitations

Expected values for the APTT test will vary from one laboratory to another, depending on the technique used. The method of clot detection, temperature, pH, collection technique, type of anticoagulant and time and method of specimen storage are all very important. Thus, laboratories should establish their own expected reference values for patients and well-defined performance standards for the control plasmas. In addition close attention should be paid to the condition of the sample. Use of hemolyzed, lipemic, or icteric samples should be avoided as these conditions may affect results, especially when using photo-optical instruments.

EXPECTED VALUES

A reference range study was conducted using frozen plasma samples from 40 normal health adults. (approximately equal numbers of males and females were used). The APTT results were as follows:
PHOSPHOLIN ES
(N= 40)

	Mean (seconds)	Range for +/- 2 SD
PHOTO-OPTICAL	29.8	23.2 – 36.4
MECHANICAL	31.8	27.1—36.5

These values should serve only as guidelines. Because differences may exist among instruments, laboratories and local populations, it is recommended that each laboratory establish its own reference range of expected APTT values.

PERFORMANCE CHARACTERISTICS

1. Precision Studies

Precision studies were performed to establish Within-Run and Between-Run CV's for normal control plasma and abnormal control plasma. A single lot number of **Phospholin ES** kits was used for these studies. Results are shown below:

Within Run (APTT in seconds)

Control	N	Mean (secs.)	S.D.	C.V.
PlasmaCon N	30	29.3	0.2	0.9%
PlasmaCon L-2	30	88.2	0.6	0.6%

Between Run (APTT in seconds)

Control	N	Mean (secs.)	S.D.	C.V.
PlasmaCon N	15	30.6	0.68	2.24%
PlasmaCon L-2	15	91.4	2.61	2.85%

II. Comparison Studies

A. Patient Testing - Comparison of **PHOSPHOLIN ES** Reagent and a commercially available reagent was performed using samples with normal and abnormal clotting activities due to either heparin therapy, factor deficiencies, or both. A total of 140 samples, measured in triplicate, were tested.

Parameters of the linear regression equations were:

Photo-optical	slope = 0.864	r = 0. 92
Mechanical	slope = 1.0929	r = 0. 93

REFERENCES:

- Proctor R.R, Rapaport S.I, The Partial Thromboplastin Time with Kaolin, AM.J.Clin. Path. 36:212-219, 1961.
- Brandt J.T., Triplett D.A., Laboratory Monitoring of Heparin. Effect of Reagents and Instruments on the APTT. AM.J.Clin. Path. 76: 530-537, 1981.
- Kelsey P.R., The Diagnosis of Lupus Anticoagulants by the Partial Thromboplastin Time, Thromb. Haemost. 52: 172-175, 1984.

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LITERATURHINWEISE:

1. Proctor R.R, Rapaport S.L., The Partial Thromboplastin Time with Kaolin. AM.J.Clin. Path. 36:212-219, 1961.
2. Brandt J.T., Triplett D.A., Laboratory Monitoring of Heparin. Effect of Reagents and Instruments on the APTT. AM.J.Clin. Path. 76: 530-537, 1981.
3. Kelsey P.R., The Diagnosis of Lupus Anticoagulants by the Partial Thromboplastin Time, Thromb. Haemost. 52: 172-175, 1984.

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	R2 Diagnostics, Inc. South Bend, Indiana USA (574) 288-4377
---	--

	Enzyme Research Laboratories, Ltd. Swansea, UK
---	---

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REFERENCIAS:

- Proctor R.R, Rapaport S.I., The Partial Thromboplastin Time with Kaolin, AM.J.Clin. Path. 36:212-219, 1961.
- Brandt J.T., Triplett D.A., Laboratory Monitoring of Heparin. Effect of Reagents and Instruments on the APTT. AM.J.Clin. Path. 76: 530-537, 1981.
- Kelsey P.R., The Diagnosis of Lupus Anticoagulants by the Partial Thromboplastin Time, Thromb. Haemost. 52: 172-175, 1984.

R2 Diagnostics, Inc.
South Bend, Indiana USA (574) 288-4377

Enzyme Research Laboratories, Ltd.
Swansea, UK

LL-4501 Rev. C