



Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK

Tel: +44 (0)1792 686070

Email: info@enzymeresearch.co.uk

Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No:	BPCa
Lot:	SAMPLE
Description:	Bovine Activated Protein C
Concentration:	0.92 mg/mL
Aliquots:	1 x 0.10 mg
Volume/aliquot:	1.087 mL
Activity:	N/A
Total Protein:	0.10 mg (1 x 0.1 mg)
Buffer:	20 mM Tris-HCl / 0.1 M NaCl / pH 7.4
Storage:	<-60 °C
Expiration Date:	** / ****

The Activated Bovine Protein C was activated from homogeneous protein C by activation with purified alpha-thrombin. The alpha-thrombin is removed after activation. Complete activation is observed on SDS-PAGE gels.

All bovine material is from U.S. cows < / = 36 months of age.

For Research Use Only.