



Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK
Tel: +44 (0)1792 686070
Email: info@enzymeresearch.co.uk
Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No:	BFIXa 1080
Lot:	SAMPLE
Description:	Bovine Factor IXa beta
Concentration:	1.00 mg/mL
Aliquots:	1 x 0.50 mg
Volume/aliquot:	0.500 mL
Activity:	3000.00 PEU/mg
Total Protein:	0.50 mg (1 x 0.5 mg)
Buffer:	20 mM Tris-HCl / 0.1 M NaCl / pH 7.4
Storage:	<-60 °C
Expiration Date:	** / ****

The Bovine Factor IXa beta was prepared from homogeneous Bovine Factor IX by activation with Bovine Factor XIa. This Bovine Factor XIa was removed after activation. Complete activation was observed on 10 % SDS-PAGE gels.

All bovine material is from U.S. cows <= 30 months of age.

For Research Use Only.