



Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK

Tel: +44 (0)1792 686070

Email: info@enzymeresearch.co.uk

Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No:	BFXIa 11a
Lot:	SAMPLE
Description:	Bovine Factor XIa
Concentration:	1.16 mg/mL
Aliquots:	1 x 1.00 mg
Volume/aliquot:	0.862 mL
Activity:	n/a
Total Protein:	1.00 mg
Buffer:	50mM Tris-HCl / 0.1 M sodium chloride / pH 8.0
Storage:	<-60 °C
Expiration Date:	**/****

The Bovine XIa is a doublet on SDS-PAGE gels and shows complete reduction upon incubation with 2-mercaptoethanol.

All bovine material is purified from U.S. cows < 30 months of age.

For Research Use Only.