



Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK
Tel: +44 (0)1792 686070
Email: info@enzymeresearch.co.uk
Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No:	BFIB 1
Lot:	SAMPLE
Description:	Bovine Fibrinogen, Plasminogen Depleted
Concentration:	45.36 mg/mL (before lyophilization)
Aliquots:	1 x 1.00 gm
Volume/aliquot:	22.046 mL
Activity:	100% Clottable
Total Protein:	1.00 gm (1 x 1.00 gm)
Buffer:	20 mM Sodium Citrate-HCl / pH 7.4
Storage:	2-8 °C
Expiration Date:	**/****

The Bovine Fibrinogen, Plasminogen depleted is >95 % pure by SDS-PAGE and is greater than 95 % clottable. The Plasminogen was depleted using a lysine-sepharose column.

All bovine material is from U.S. cows <=30 months of age.

Hydrate the protein to the original volume with warmed distilled water/buffer. The hydration should take place in a 37 °C water bath to ensure all protein solubilizes. After hydration, aliquot and freeze the unused portions at <-60 °C.

For Research Use Only.