

Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK
 Tel: +44 (0)1792 686070
 Email: info@enzymeresearch.co.uk
 Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No: FVII-DP
Lot: SAMPLE
Product: Factor VII Deficient Human plasma.
Description: Pooled normal citrated human plasma depleted of FVII using antibodies directed to FVII immobilized on agarose beads. Plasma contains 20 mM HEPES.
Presentation: Frozen Factor VII deficient plasma.
Preparation/Handling: Thaw 1 mL vials in 37°C waterbath for 5 minutes; for bulk volumes, thawing time will be dependant on bottle size.
Storage and Stability: Plasma is shipped frozen and should be stored below -60°C. Product is stable until date stated on vial label when stored at -60°C. Once thawed, plasma is stable for 4 hours at 2-8°C in original vial.
Expiration Date: **/****
Specifications: Lot specific assay results;

Activity Assay	
PT (Stago Neoplastine C1+)	45.6 sec
APTT (Hemosil APTT-SP)	31.3 sec
Fibrinogen (clottable)	3.25 g/L
FII (activity)	1.08 U/mL
FV (activity)	0.95 U/mL
FVII (activity)	<0.01U/mL
FVIII (activity)	0.98 U/mL
FIX (activity)	1.39 U/mL
FX (activity)	1.09 U/mL
FXI (activity)	1.29 U/mL
FXII (activity)	1.46 U/mL
FVII Inhibitor assay (1:1, 120 min @37°C)	negative

Although this material was prepared with plasma collected from donors screened for CJD and which was tested at source and found negative for HBsAg, Syphilis and antibodies to HIV and HCV and non-reactive for HIV-1 rNA and HCV rNA by FDA approved tests, it should be handled by personnel trained in the proper procedures for handling potential viral contaminants.

Related Products

SAFVII-IG	Sheep Anti-Human FVII, whole IgG from antiserum	10 mg vial
SAFVII-AP	Sheep Anti-Human FVII, affinity purified IgG	0.5 mg vial
SAFVII-HRP	Sheep Anti-Human FVII, AP IgG-peroxidase conjugate	0.2 mg vial
SAFVII-APFTC	Sheep Anti-Human FVII, FITC-conjugated affinity purified IgG	0.1 mg vial
SAFVII-APBIO	Sheep Anti-Bovine FVII, biotinylated affinity purified IgG	0.1 mg vial
FVII-EIA	Paired Antibody set for ELISA of FVII antigen	5 x 96 wells

For Research Purposes Only