



Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK  
Tel: +44 (0)1792 686070  
Email: [info@enzymeresearch.co.uk](mailto:info@enzymeresearch.co.uk)  
Website: [www.enzymeresearch.co.uk](http://www.enzymeresearch.co.uk)

## CERTIFICATE OF ANALYSIS

Cat No:	<b>Mab HFIX</b>
Lot:	SAMPLE
Description:	Monoclonal antibody to Human Factor IX
Concentration:	1 mg/mL
Aliquots:	1 x 0.50 mg
Volume/aliquot:	0.500 ml
Activity:	N/A
Total Protein:	0.50 mg
Buffer:	10 mM Sodium Phosphate / 0.15 M NaCl / 0.01% Sodium Azide / pH 7.4
Storage:	-20°C
Expiration Date:	**/****

Epitope on Factor IX Light Chain.

Antibody purified by Protein G affinity chromatography from cell supernatant. Binds human factor IX in solid phase ELISA. Blots human factor IX light chain.

We suggest hydrating the antibody to its original concentration of 1 mg/mL and aliquot to single use size (avoid freeze/thaw cycles). Antibody expires 2 years after reconstitution.

***For Research Use Only.***