



Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK

Tel: +44 (0)1792 686070

Email: info@enzymeresearch.co.uk

Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No:	RFXa
Lot:	SAMPLE
Description:	Rat Factor Xa
Concentration:	2.14 mg/mL (before lyophilization)
Aliquots:	1 x 0.10 mg
Volume/aliquot:	0.047 mL
Activity:	N/A
Total Protein:	0.10 mg
Buffer:	20 mM Tris-HCl, 0.7 M NaCl, pH 7.4
Storage:	2-8 °C
Expiration Date:	**/****

Rat Factor Xa is prepared from homogeneous Rat Factor X by activation with Russells' Viper Venom. This RVV-X is removed after activation. Complete activation is observed on SDS-PAGE. We suggest hydrating the protein to the original volume with water.

For Research Use Only.