Enzyme Research Laboratories
Global Coag

Glamorgan House . Charter Court . Enterprise Park . Swansea . SA7 9DD . UK
Tel: +44 (0)1792 686070
Email: info@enzymeresearch.co.uk
Website: www.enzymeresearch.co.uk

CERTIFICATE OF ANALYSIS

Cat No:	**BFXa**
Lot:	SAMPLE
Description:	Bovine Factor Xa
Concentration:	1.01 mg/mL
Aliquots:	1 x 1.00 mg
Volume/aliquot:	0.990 mL
Activity:	249.78 IU/mg
Total Protein:	1.00 mg (1 x 1.0 mg)
Buffer:	20 mM Tris-HCl / 0.7 M NaCl, pH 7.4
Storage:	<-60°C
Expiration Date:	**/****

The Bovine Factor Xa was prepared from homogeneous Bovine Factor X by activation with Russells' Viper Venom. This RVV-X is removed after activation. Complete activation was observed on SDS-PAGE. Product is >/=95% pure upon visual inspection of SDS-PAGE.

All bovine material is from USDA inspected herd pre and post mortem. All Bovine material is from U.S. cows </=30 months of age.

For Research Use Only.